DESCRIPTION

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Instructions for Use

DESCRIPTION

i=FACTOR™ Bone Graft is a pure, natural form of hydroxyapatite (anorganic bone mineral or ABM) containing P-15™, a synthetic peptide. It is provided in two forms, Putty and Flex. Putty is comprised of ABM/P-15 particles that are suspended in an inert biocompatible hydrogel. Flex is a lyophilized (freeze-dried) version of Putty that is formed into strips.

Hydroxyapatite is the major mineral component of bone. P-15 is a fifteen amino acid sequence polypeptide, a portion of Type-I collagen that serves as a biomimetic attachment site for bone cells. The hydrogel carrier (glycerin and carboxymethylcellulose or CMC) forms the matrix that improves the overall handling characteristics and helps contain the i=FACTOR Bone Graft at the surgical site. The bovine-derived ABM particles within the matrix are radiopaque and sized between 250 and 425 microns. i=FACTOR Bone Graft provides a bone graft substitute that is remodeled into new bone during the natural healing process.

INDICATIONS FOR USE

i=FACTOR Bone Graft is intended to fill bony voids or gaps in extremities and spine that are not intrinsic to the stability of the bony structure. Sufficient internal or external fixation is required.

INTENDED USER

i=FACTOR Bone Graft is a bone graft substitute material for use only by a qualified Orthopedic surgeon, Neurosurgeon, or Spine surgeon.

CONTRAINDICATIONS

Contraindications include but are not limited to:

- Absence of load bearing structural support at the graft site
- Sensitivity to components of the i=FACTOR Bone Graft
- Active infection at the operative site
- Operative site subject to excessive impact or stress
- Significant vascular impairment proximal to the graft site
- Use in direct contact with articular spaces
- Presence of segmental defects
- Metabolic or systemic bone disorders that affect bone or wound healing
- Compromised renal function
- Uncooperative patients who will not or cannot follow postoperative Instructions, including individuals who abuse drugs and/or alcohol

WARNINGS AND PRECAUTIONS

i=FACTOR Bone Graft is not intended to provide load bearing structural support during the healing process. Rigid fixation techniques are recommended as needed to assure stabilization of the defect in all planes. As with any surgical procedure care should be exercised in treating individuals with pre-existing conditions that may affect the success of the surgical procedure. This includes, but is not limited to, individuals with bleeding disorders of any etiology; calcium metabolism disorders, e.g., hyperglycemia; long-term steroidal therapy; immunosuppressive therapy; or high dosage radiation therapy.

DO NOT USE IF STERILE PACKAGING IS OPENED OR DAMAGE. Discard or return damaged packaging and all contents.

i=FACTOR Bone Graft is designed for single patient use only. Do not attempt to re-sterilize or re-use. Discard unused contents. Attempting to reuse i=FACTOR Bone Graft will adversely affect product sterility and physical handling characteristics.

i=FACTOR Bone Graft should only be used in surgical procedures where it can be adequately contained at the bony void or defect. Avoid overfilling the bone void or pressurizing the treatment site. Inadequate containment of i=FACTOR Bone Graft could result in product migration from the intended bony defect site. If product migration occurs, clinical outcomes may be compromised by the lack of bone graft material in the appropriate space. Potential patient adverse events caused by inadequate containment and migration of i=FACTOR Bone Graft could include, but are not limited to the following: pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function; any of which may require revision surgery.

The effect of i=FACTOR Bone Graft on pregnant or nursing patients has not been evaluated.

The use of i=FACTOR Bone Graft when mixed with other bone graft substitute products has not been evaluated; therefore, the effectiveness of i=FACTOR when used in this manner is unknown.

ADVERSE EFFECTS

Possible adverse effects include, but are not limited to:

- Wound complications including hematoma, site drainage, infection and other complications that are possible with any surgery
- Extrusion or migration of the bone void filler, as is possible with any bone void filler, resulting in pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function; any of which may require revision surgery
- Nonunion, malunion or delayed union
- Loss of reduction
- Refracture
- Incomplete, or lack of osseous ingrowth into the bone void, as is possible with any bone void filler
- Transient hypercalcemia
- Allergic reaction to components of the i=FACTOR Bone Graft

STORAGE

i=FACTOR Bone Graft in its original packaging can be stored at ambient room temperature.

DIRECTIONS FOR USE

PUTTY

Remove syringe cap and dispense i=FACTOR Bone Graft by depressing the syringe plunger. i=FACTOR Bone Graft may be dispensed directly to the surgical site or into a separate sterile receptacle where it can be transferred using traditional surgical instrumentation or by hand.

i=FACTOR Bone Graft should only be placed in voids where it can be contained adequately.

NOTE: When opening the foil pouch containing the i=FACTOR Bone Graft syringe, you may notice a very small amount of water retained within the pouch. This is a normal part of the steam sterilization process and does not affect the integrity or sterility of the product.

FLEX

Remove i=FACTOR Bone Graft from the tray. i=FACTOR Bone Graft may be placed directly to the surgical site or into a separate sterile receptacle where it can be transferred using traditional surgical instrumentation or by hand.

i=FACTOR Bone Graft should only be placed in voids where it can be contained adequately.

 MANUFACTURED BY:

Cerapedics, Inc.
11025 Dover Street, Suite 1600
Westminster, CO 80021 USA
+1 303.974.6275 Voice
+1 303.974.6285 Fax