Instructions for Use

Bone Graft (i-FACTOR)

Disposable Sterilized Medical Device (Do Not Reuse)

- Product Name: i-FACTOR Putty
- Item Name: Bone Graft
- Model Name: Refer to the label on the side of this box [REF]
- Expiration Date: Refer to the label on the side of this box [☞ Year-Month-Date]
- Importer: SpinoTech (5th F/L, Sehyun Bldg 7, Saimdang-ro-1-gil, Seocho-gu, Seoul)
- Manufacturer: Cerapedics, Inc. / Company: U.S.A
- Product License Number: Import 18-362
- Manufacturing number: Refer to the label on the side of this box [LOT]
- Packing: 1 Unit
- Purpose of Use: This product is used for spinal fusion. It should be used in combination with supplemental plate fixation in skeletally mature patients.
- Storage condition: Store at dry room temperature
- Sterilization method: Steam Sterilization

■ **How to use**

A. Pre-Use Ready

- 1) Check the expiration date, appearance, and packaging of the product before use.
- 2) Read the manufacturer's IFU before using it.

B. How to use

- 1) Remove the syringe cap and push the plunger to inject the i-FACTOR Bone Graft.
 - 2) To inject the i-FACTOR Bone Graft directly into the surgical site or with a normal surgical instrument or hand.
 - 3) Use an appropriate amount for the bone loss area.
 - 4) This product must be used in combination with supplemental plate fixation such as intervertebral fixation and metal bone fixation licensed by the KFDA. The position is evaluated during surgery with simple radiographic or fluoroscopy.
- ※ When opening the foil pouch containing the i-FACTOR Bone Graft syringe, you may notice a very small amount of water retained within the pouch. This is a normal part of the steam sterilization process and does not affect the integrity or sterility of the product.
- ※ The use of this product for patients under the age of 18 can have a negative effect on bone development by forming excessive bone and changing growth. Therefore, use is prohibited for patients under the age of 18 with severe bone defects or rare bone disorders.

C. How to store and manage after use

This product is a disposable sterilized product and is prohibited from being reused.

■ **Precaution for use**

A. Contraindications

- 1) Absence of load bearing structural support at the graft site
- 2) Sensitivity to components of this product

- 3) Active infection at the operative site
- 4) Operative site subject to excessive impact or stress
- 5) Significant vascular impairment proximal to the graft site
- 6) Use in direct contact with articular spaces
- 7) Presence of segmental defects
- 8) Metabolic or systemic bone disorders that affect bone or wound healing
- 9) Compromised renal function
- 10) Uncooperative patients who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol
- 11) Immune system damage

B. Warnings AND Precautions

- 1) This product is not intended to provide load bearing structural support during the healing process. Rigid fixation techniques are recommended as needed to assure stabilization of the defect in all planes. As with any surgical procedure, care should be exercised in treating individuals with pre-existing conditions that may affect the success of the surgical procedure. This includes, but is not limited to, individuals with bleeding disorders of any etiology, long term steroidal therapy, immunosuppressive therapy or high dosage radiation therapy.
- 2) Do not use if sterile packaging is opened or damaged. Discard or return damaged packaging and all contents.
- 3) This product is designed for single patient use only. Do not attempt to re-sterilize or re-use. Discard unused contents. Attempting to reuse this product will adversely affect product sterility and physical handling characteristics.
- 4) This product should only be used in surgical procedures where it can be adequately contained at the bony void or defect. Avoid overfilling the bone void or pressurizing the treatment site.
- 5) Inadequate containment of this product could result in product migration from the intended bony defect site. If product migration occurs, clinical outcomes may be compromised by the lack of bone graft material in the appropriate space.
- 6) Potential patient adverse events caused by inadequate containment and migration of this product could include, but are not limited to the following: pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function; any of which may require revision surgery.
- 7) The effect on pregnant or nursing patients has not been evaluated.
- 8) The use of this product when mixed with other bone graft substitute products has not been evaluated. Therefore, the effectiveness of this product when used in this manner is unknown.
- 9) The use of this product for patients under the age of 18 can have a negative effect on bone development by forming excessive bone and changing growth. Therefore, use is prohibited for patients under the age of 18 with severe bone defects or rare bone disorders.
- 10) Safety in application areas where strength is required has not been confirmed.
- 11) This product must be used in combination with supplemental plate fixation such as intervertebral fixation and metal bone fixation licensed by the KFDA. The position is evaluated during surgery with simple radiographic or fluoroscopy.

C. Adverse Effect

- 1) Transient hypercalcemia
- 2) Allergic reaction to components of this product

D. Adverse Event

- 1) Wound complications including hematoma, site drainage, infection and other complications that are possible with any surgery
- 2) Nonunion, malunion or delayed union
- 3) Incomplete or lack of osseous ingrowth into the bone void, as is possible with any bone void filler
- 4) Extrusion or migration of the bone void filler, as is possible with any bone void filler, resulting in pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function; any of which may require revision surgery
- 5) Refracture
- 6) Loss of reduction

■ **Related Institution on Adverse Event Reporting**

NIDS(National Institute of Medical Device Safety Information), 080-080-4183

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