iFACTOR+ Matrix Bone Graft Instructions for Use

DESCRIPTION

iFACTOR+ Matrix is a bone graft substitute material comprised of porcine-derived anorganic bone mineral particles (pABM) coated with P-15, a synthetic peptide. The pABM/P-15 particles are embedded in a carrier consisting of bovine-derived Type I dermal collagen that aids in the handling and containment of the particles within the intended surgical site. Upon implantation in a bony site, iFACTOR+ Matrix is remodeled into new native bone by cell-mediated resorption during the natural healing process.

The pABM component of iFACTOR+ Matrix is a natural form of hydroxyapatite (calcium phosphate), which is the main mineral component of bone. P-15 is a synthetically derived fifteen amino acid polypeptide that mimics the cell binding domain of human Type I collagen, and serves as an attachment site for osteogenic cells. The pABM/P-15 particles are naturally porous, radiopaque, and sized at 250-1000 microns (nominal average = 481 microns). iFACTOR+ Matrix is composed of 80% pABM/P-15 particles and 20% bovine collagen.

iFACTOR+ Matrix is provided as a freeze-dried material that, when hydrated, forms a moldable putty that can be shaped as desired. iFACTOR+ Matrix is provided sterile and is intended for single use only.

INDICATIONS FOR USE

iFACTOR+ Matrix Bone Graft is intended to fill bony voids or gaps in extremities and spine that are not intrinsic to the stability of the bony structure. Sufficient internal or external fixation is required.

iFACTOR+ Matrix may be mixed with autogenous bone and/or autogenous bone marrow aspirate (BMA).

INTENDED USER

iFACTOR+ Matrix is a bone graft substitute material for use only by a qualified Orthopedic surgeon, Neurosurgeon, or Spine surgeon.

CONTRAINDICATIONS

Use of iFACTOR+ Matrix is contraindicated for any of the following:

- Absence of load bearing structural support at the graft site
- Patients with known sensitivity to any components of iFACTOR+ Matrix
- Active infection at the operative site
- Operative site subject to excessive impact or stress
- Significant vascular impairment proximal to the graft site

- Use in direct contact with articular spaces
- Presence of segmental defects
- Metabolic or systemic disorders that affect bone or wound healing
- Compromised renal function
- Uncooperative patients who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol

WARNINGS AND PRECAUTIONS

- i≡FACTOR+ Matrix is not intended to provide load-bearing structural support during the healing process. Rigid fixation techniques are recommended as needed to ensure stabilization of the bony structure in all planes.
- i≡FACTOR+ Matrix should only be used in surgical procedures where it can be adequately contained at the bony void or defect. Inadequate containment of i≡FACTOR+ Matrix could result in product migration from the intended bony defect site. If product migration occurs, clinical outcomes may be compromised by the lack of bone graft material in the appropriate space. Potential patient adverse events caused by inadequate containment and migration of i≡FACTOR+ Matrix may include, but are not limited to the following: pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function; any of which may require revision surgery.
- Avoid overfilling the bone void or pressurizing the treatment site.
- As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure. This includes, but is not limited to, individuals with bleeding disorders of any etiology, long-term steroidal therapy, immunosuppressive therapy or high dosage radiation therapy.
- The use of i≡FACTOR+ Matrix when mixed with other bone graft substitute products has not been evaluated; therefore, the effectiveness of i≡FACTOR+ Matrix when used in this manner is unknown.
- The effect of i≡FACTOR+ Matrix on pregnant or nursing patients has not been evaluated.
- **DO NOT USE IF STERILE PACKAGING IS OPENED OR DAMAGED.** Discard or return damaged packaging and all contents.
- i≡FACTOR+ Matrix is designed for single patient use only. Do not attempt to re-sterilize or re-use.
- Do not use after the printed expiration date on the label.
- Discard any unused graft material.

POTENTIAL ADVERSE EVENTS

As with any bone graft material, the following potential adverse events may occur. Occurrence of one or more of these conditions may require an additional surgical procedure and may also require removal of i≡FACTOR+ Matrix.

- Extrusion or migration of the bone graft material resulting in pain, neural impingement, physical impairment, irritation or wear of an articulating joint, or loss of function
- Wound complications including hematoma, superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, cellulitis, site drainage, or wound dehiscence
- Incomplete or lack of osseous ingrowth into the graft site
- Excessive or abnormal bone formation in an unintended location
- Delayed union, nonunion, or malunion
- Loss of reduction, refracture, cyst recurrence
- Arthritis or other disorders in bone formation
- Transient hypercalcemia
- Allergic reaction to components of i≡FACTOR+ Matrix
- Localized immunological reactions consisting of transient localized edema, swelling, and rash have been reported to occur with bone void fillers containing collagen. Although there is no evidence that the device will be unsafe or ineffective in such patients, the safety and effectiveness of the device in these patients has not been established.
- Renal problems
- Nervous system problems

STORAGE

Store i≡FACTOR+ Matrix in its original packaging at ambient room temperature. Do not freeze or expose to extreme heat.

DIRECTIONS FOR USE

Familiarization with i≡FACTOR+ Matrix and proper bone grafting and rigid fixation techniques are extremely important. Radiographic evaluation of the planned graft site is essential for accurately assessing the desired surgical outcome, and to aid in the selection and placement of i≡FACTOR+ Matrix along with any fixation devices.

To prepare i≡FACTOR+ Matrix, hydrate the i≡FACTOR+ Matrix material with either sterile saline or autogenous bone marrow aspirate (BMA) using the following recommended fluid volumes:

i≡FACTOR+ Matrix	Fluid volume (saline or BMA)
1.0 cc	1.0 cc
2.5 cc	2.4 cc
5.0 cc	4.4 cc
10.0 cc	9.6 cc

Once i=FACTOR+ Matrix has hydrated, knead the material to a putty-like consistency for at least 60 seconds and mold into the desired shape. Insert the bone graft into the surgical site.

Fixation of the implant site must be sufficient to prevent collapse and deformity secondary to functional loading. Anatomical reduction and rigid fixation in all planes must be obtained to ensure that the graft is not supporting load.

Postoperative patient management should follow the same regimen as with other bone grafts or autogenous bone graft. Standard postoperative practices should be followed, particularly as applicable to surgical procedures involving the use of fixation devices.

WARRANTIES

All warranty rights are lost if modifications are made to this product. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

PRODUCT COMPLAINTS

Any health care professional using i=FACTOR+ Matrix who has complaints about the product, or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Cerapedics. Further, if any of the implanted product ever “malfunctions”, i.e., does not meet any of its performance specifications or otherwise does not perform as intended, or may have caused or contributed to the death or serious injury of a patient, Cerapedics should be notified immediately by telephone, or through written correspondence via email or fax. When filing a complaint, please provide the product name, the product reference code, lot number, your name and address, and the nature and details of the complaint.

FURTHER INFORMATION

If further information is required, please contact Cerapedics at the address below.

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